

Aesculap Orthopaedics

Columbus[®] MIOS Tibia Plateau

Minimally Invasive
Orthopaedic Solutions



Product Information / Surgical Technique Supplement

The minimally invasive tibia plateau, with proven stability




The Columbus® MIOS Tibial Plateau is designed to facilitate implantation through the smaller incisions used for minimally invasive knee replacement surgery. Together with the OrthoPilot® Computer-Assisted Navigation System, the Columbus® MIOS System provides one of the most comprehensive solutions to minimally invasive techniques available today.

The special design of the Columbus® MIOS Tibial Plateau offers the following advantages:

- Easy and safe implantation due to the small and optimized stem
- Minimized potential for malalignment
- Helps avoid damage and trauma to the soft tissue due to overstretching
- Reduction of the risk of femoral condyle damage
- The specially designed MIOS instruments helps lead to a safe procedure for the patient
- MIOS approach preserves bone for the future


The Columbus® MIOS System can lead to an improvement of the functional outcomes (especially in the early postoperative period), reduce the perioperative morbidity, and accelerate postoperative recovery.

Let us break new ground together.



Proximal Distal Micro Motion




- Test performed at 2.7 kN axial load (average axial load while walking)
- Proximal distal micro motion while walking (maximum tolerated motion: 150 μm)
- Worst case simulation tests plateaus with 15° anterior slope
- Outcome based on 3,000 cycles



Surgical Technique Supplement

How it works

- Position the proper trial tibia plateau on the resected tibia in the correct ML and AP position. Insert the pre-fixation peg assembled with the holder (fig. 1).
- Perform trial after removing the holder and pre-fixation guide to establish the natural rotational position under the femoral trial component. Fix the selected position with two short threaded fixation pins (fig. 2).
- Remove the pre-fixation peg. Prepare the proximal tibia with the appropriate size broach (fig. 3).



Columbus® MIOS Tibia Plateau CR/PS –ST0255

Uncoated


NN370K	Columbus® CR/PS MIOS Tibia Plateau T0
NN368K	Columbus® CR/PS MIOS Tibia Plateau T0+
NN371K	Columbus® CR/PS MIOS Tibia Plateau T1
NN372K	Columbus® CR/PS MIOS Tibia Plateau T1+
NN373K	Columbus® CR/PS MIOS Tibia Plateau T2
NN374K	Columbus® CR/PS MIOS Tibia Plateau T2+
NN375K	Columbus® CR/PS MIOS Tibia Plateau T3
NN376K	Columbus® CR/PS MIOS Tibia Plateau T3+
NN377K	Columbus® CR/PS MIOS Tibia Plateau T4
NN378K	Columbus® CR/PS MIOS Tibia Plateau T4+
NN379K	Columbus® CR/PS MIOS Tibia Plateau T5



Overview of the most important dimensions for Columbus® MIOS Tibia Plateau


	T0/T0+	T1/T1+	T2/T2+	T3/T3+	T4/T4+	T5
A	62	65	70	75	80	85
B	41/44	43/46	45/49	48/52	51/55	56
C	14/14.5	15/16	16/17.5	17.5/19	19/20.5	20.5
D	12.3	12.3	12.3	12.3	14.3	14.3

measurements in [mm]





Columbus® MIOS Tibia Instrument Set—ST0254





AESCULAP[®]

Implant Systems

Aesculap Implant Systems, Inc.

3773 Corporate Parkway
Center Valley, PA 18034
Phone: 866-229-3002

www.aesculapimplantsystems.com

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2008 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA. DOC751 3M 2/08
Aesculap is an equal opportunity employer