



Aesculap Endoscopic Technology

Aesculap® AdTec® mini Single Use

3.5 mm single-use instruments for minimal trauma

3,5 mm Einmal-Instrumente für minimales Trauma

NEW

Mini Laparoscopy with AdTec® mini

- Small incision
- Less post-operative pain¹
- Less visible scars¹
- No change of laparoscopic operating technique

Mini Laparoskopie mit AdTec® mini

- Kleine Hautinzision
- Weniger postoperative Schmerzen¹
- Weniger sichtbare Narben¹
- Keine Änderung der laparoskopischen OP Technik








Statements relate to Mini Laparoscopy / Angaben beziehen sich auf Mini Laparoskopie

¹ Novitsky Y, Kercher K W, Czerniach D R, Kaban G K, Khera S, Gallagher-Dorval K A, Callery M P, Litwin D E M, Kelly J J.

Advantages of Mini-laparoscopic vs Conventional Laparoscopic Cholecystectomy. Arch Surg. 2005 Dec;140(12):1178-83.

Aesculap® AdTec® mini Single Use

		Diameter Ø	length Länge
AdTec® mini Single use <ul style="list-style-type: none"> Ergonomic handle with large contact surface for fingers One-to-one tactile feedback Shaft 360° rotatable Integrated overload protection in handle No protruding joints distally Insulated to the tips Pack containing 5 pieces 	PO830SU Metzenbaum scissors Metzenbaum Schere	3.5 mm	290 mm
			
	PO831SU Maryland grasping and dissecting forceps Maryland Fass- und Präparierzange	3.5 mm	290 mm
AdTec® mini Single use <ul style="list-style-type: none"> Ergonomischer Handgriff mit großer Auflagefläche für die Finger Hervorragende Taktilität Schaft 360° drehbar Integrierter Überlastschutz im Handgriff Keine austretende Gelenke Isolierung bis an die Maulteile Packung á 5 Stück 			
	PO833SU Atraumatic universal grasping forceps Atraumatische Fasszange	3.5 mm	290 mm
			
Aesculap® MiniPort <ul style="list-style-type: none"> Compatible single use 3.5 mm trocar Kompatibler 3,5 mm Einwegtrokar 	EK105SU MiniPort Disposable Trocar Sleeve, with integrated seal and thread diam.: 3.5 mm, length: 70 mm, Package á 6 pcs. MiniPort Einmal-Trokarhülse mit Dichtung und Gewinde, Ø 3,5 mm, Länge: 70 mm, Packung á 6 Stück		
			
	EK030R Dilating Trocar Pin, diam.: 3.5 mm Dilatations Trokardorn, Ø 3,5 mm		
			

The main product trademark 'Aesculap' is a registered trademark of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.