


Aesculap Orthopaedics

Excia[®] L


Lateralized Hip Stem



- Increased joint stability without affecting leg length
- Approximately 6mm additional lateral offset
- Increased varus neck angle of 128°
- Seamless integration with OrthoPilot[®] Navigation System

Excia® L

The Excia L stems offer additional options for tensioning the joint and controlling leg length. Combined with Aesculap's OrthoPilot® Navigation System, the surgeon can determine the exact leg length and offset for all 10 head positions during the entire rasping procedure, allowing adjustments to be made even before trial reduction of the joint. Contact Aesculap Implant Systems to learn more about Excia L.



Excia L Implants	
Item #	Description
NC428T	Plasmapore μ-CaP 8/10 8mm
NC429T	Plasmapore μ-CaP 8/10 9mm
NC430T	Plasmapore μ-CaP 8/10 10mm
NC431T	Plasmapore μ-CaP 8/10 11mm
NC432T	Plasmapore μ-CaP 8/10 12mm
NC433T	Plasmapore μ-CaP 8/10 13mm
NC434T	Plasmapore μ-CaP 8/10 14mm
NC435T	Plasmapore μ-CaP 8/10 15mm
NC436T	Plasmapore μ-CaP 8/10 16mm
NC437T	Plasmapore μ-CaP 8/10 17mm
NC438T	Plasmapore μ-CaP 8/10 18mm

*Sizes 8mm-11mm not yet available in the USA

Excia Complete Instrument Set # UQ0203			
Excia General Instrument Tray		Excia Rasp Tray	
Item #	Description	Item # Description	
NG901R	Tray for Excia Instruments	NG903R	Tray for Excia Rasps
JF217R	Tray Lid	JF217R	Tray Lid
ND004R	Femoral Head Saw Guide	ND769R	Excia Box Osteotome
ND050	Head Impactor	NG906R	Excia Wing Rasp
ND375	Head Remover	NG908R	Excia Rasp Size 8
ND401R	Implant Removal Punch	NG909R	Excia Rasp Size 9
NF340R	Stem Impactor	NG910R	Excia Rasp Size 10
NG181	Trial Head Lateralized 28mm S	NG911R	Excia Rasp Size 11
NG182	Trial Head Lateralized 28mm M	NG912R	Excia Rasp Size 12
NG183	Trial Head Lateralized 28mm L	NG913R	Excia Rasp Size 13
NG184	Trial Head Lateralized 28mm XL	NG914R	Excia Rasp Size 14
NG185	Trial Head Lateralized 28mm XXL	NG915R	Excia Rasp Size 15
NG191	Trial Head Lateralized 32mm S	NG916R	Excia Rasp Size 16
NG192	Trial Head Lateralized 32mm M	NG917R	Excia Rasp Size 17
		NG918R	Excia Rasp Size 18
		NG929R	Excia Rasp Handle for Lat. Pos.
NG193	Trial Head Lateralized 32mm L		
NG194	Trial Head Lateralized 32mm XL		
NG195	Trial Head Lateralized 32mm XXL		
NG301	Trial Head 28mm S		
NG302	Trial Head 28mm M		
NG303	Trial Head 28mm L		
NG304	Trial Head 28mm XL		
NG305	Trial Head 28mm XXL		
NG316	Trial Head 32mm S		
NG317	Trial Head 32mm M		
NG318	Trial Head 32mm L		
NG319	Trial Head 32mm XL		
NG320	Trial Head 32mm XXL		

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.

AESCULAP
Implant Systems

Aesculap, Inc.

3773 Corporate Parkway
Center Valley, PA 18034
Phone: 800-282-9000

<http://www.aesculapusa.com>